File System Command Group

File System Overview

You use the commands in the File System Command Group to help you use the built-in hard disk drive and floppy disk drive. You can use the commands to do the following:

- List the contents of the default directory
- Create and delete directories
- Create, copy, read, rename, or delete a file

File System Commands

Command	Description
FILESystem?	Returns the file system state
FILESystem:COPy	Copies one or more files to a new file
FILESystem:DELEte	Deletes a named file or directory
FILESystem:DIR?	Returns a list of directory contents
FILESystem:MKDir	Makes a new directory
FILESystem:PRInt	Prints a named file to the named port
FILESystem:READFile	Copies the named file to the GPIB port
FILESystem:REName	Assigns a new name to an existing file
FILESystem:RMDir	Deletes the named directory
FILESystem:WRITEFile	Copies the GPIB port block data to a named file

FILESystem?

Description

This query-only command returns the state of the file system.

Group

File System

Related Commands

FILESystem:COPy (see page 135), FILESystem:DELEte (see page 136), FILESystem:DIR? (see page 137), FILESystem:PRInt (see page 139), FILESystem:READFile (see page 140), FILESystem:REName (see page 141), FILESystem:WRITEFile (see page 143)

Syntax

FILESystem?

Arguments

None.

Example

FILESystem?

This query might return : FILESYSTEM: DIR, indicating that a list of files and directories exist in the default directory.

FILESystem:COPy

Description

This command (no query form) copies a named file or files to a new file. The new file may be in a totally separate directory than the old file. Also, you can use wild card characters (*.*) to copy multiple files with one command.

Group

File System

Related Commands

FILESystem: DELEte (see page 136)

Syntax

```
FILESystem:COPy {<source file path>,
<destination file path>|<source directory path>,
<destination file path>|<source directory path>,
<destination directory path>}
```

Arguments

• <file path>

This is a quoted string that defines the file name and path. Input the file path using the form <drive>/<dir>/<filename>.<drive> and one or more <dir>s are optional. <filename> stands for a filename of up to 128 characters. While filename extensions are not required, they are highly recommended. You can also use the inclusive filename *.* in the source file path to copy all files.

• <directory path>

This is a quoted string that defines the directory. Input the directory using the form <drive>/<dir></directory name>.<drive> and one or more <dir>s are optional.</directory name> stands for a directory name of up to 128 characters. While filename extensions are not required, they are highly recommended.

Example

```
FILESystem: COPy C:/MYDIR/"TEK00001,SET", "C:/TEK00001.SET"
```

This command copies the file named TEK00001.SET, located in the MYDIR directory on the C drive to a file named TEK00001.SET in the root directory.

FILESystem: DELEte

Description

This command (no query form) deletes a named file or directory.

Group

File System

Related Commands

FILESystem:COPy (see page 135)

Syntax

FILESystem: DELEte <file path>

Arguments

• <file path>

This is a quoted string that defines the file name and path. Input the file path using the form <drive>/<dir>/<filename>.<drive> and one or more <math><dir>>s are optional. If you do not specify them, the oscilloscope will delete the file in the current directory. <filename> stands for a filename of up to 128 characters. While filename extensions are not required, they are highly recommended. You can also use the inclusive filename *.* to delete all files.

Example 1

FILESystem: DELEte "NOT_MINE.SET"

This command deletes the file named NOT MINE.SET on the default drive and directory.

FILESystem:DIR?

Description

This query-only command returns a list quoted string. Each string contains the name of a file or directory in the current directory.

Group

File System

Related Commands

FILESystem: MKDir (see page 138)

Syntax

FILESystem:DIR?

Arguments

None.

Example

FILESystem:DIR?

This query returns a list of files and directories in the default directory.

FILESystem: MKDir

Description

This command (no query form) creates a new directory.

Group

File System

Related Commands

FILESystem:DIR? (see page 137)

Syntax

FILESystem: MKDir <directory path>

Arguments

<directory path>

This is a quoted string that defines the directory. Input the directory using the form <drive>/<dir></directory name>.<drive> and one or more <dir>s are optional.</directory name> stands for a directory name of up to 128 characters. While filename extensions are not required, they are highly recommended.

The current directory refers to the name of a directory as returned to the FILESystem:CWD? query.

Example

FILESystem: MKDir C:/"NEATPICS"

This command creates the directory named ${\tt NEATPICS}$ on the current drive.

FILESystem:PRInt

Description

This command (no query form) prints a named file to a named port.

Group

File System

Related Commands

FILESystem: READFile (see page 140)

Svntax

FILESystem: PRInt <filepath>, GPIb

Arguments

• <file path>

This is a quoted string that defines the file name and path. Input the file path using the form <drive>/<dir>/<filename>.<drive> and one or more <math><dir>> are optional. <filename> stands for a filename of up to 128 characters. While filename extensions are not required, they are highly recommended. You can also use the inclusive filename *.* to print all files.

• GPIh

This specifies that the hardcopy is sent out the GPIB port.

Example

FILESystem:PRInt C:/"TEK00000.IBM",GPIb

This command sends the file named $\mathtt{TEK00000}$. IBM out the GPIB port.

FILESystem:READFile

Description

This command (no query form) copies a named file to the GPIB port.

Group

File System

Related Commands

FILESystem:PRInt (see page 139)

Syntax

FILESystem: READFile <filepath>

Arguments

• <file path>

This is a quoted string that defines the file name and path. Input the file path using the form <drive>/<dir>/<filename>.<drive> and one or more <dir>s are optional. <filename> stands for a filename of up to 128 characters. While filename extensions are not required, they are highly recommended. You can also use the inclusive filename *.* to read all files.

Example

FILESystem:READFile C:/"ABC.DEF"

This command reads a file named "ABC.DEF" on the C drive and sends it to the GPIB port.

FILESystem:REName

Description

This command (no query form) assigns a new name to an existing file.

Group

File System

Syntax

FILESystem: REName <old file path>, <new file path>

Arguments

• <old file path>

This is a quoted string that defines the file to rename. Input the file path using the form <drive>/<dir>/<filename>.<drive> and one or more <dir>s are optional. <filename> stands for a filename of up to 128 characters. While filename extensions are not required, they are highly recommended.

<new file path>

This is a quoted string that defines the new name of the file. Input the file path using the form <drive>/<dir>/<filename>.<drive> and one or more <dir>s are optional. <filename> stands for a filename of up to 128 characters. While filename extensions are not required, they are highly recommended.

Example

FILESystem: REName C:/"TEK00000.SET", C:/"MYSETTING.SET"

This command gives the file named TEK000000. SET the new name of MYSETTING. SET. The file remains in the root directory on the C drive.

FILESystem:RMDir

Description

This command (no query form) deletes a named directory. The directory and all of the files and subdirectories it contains are deleted.

Group

File System

Syntax

FILESystem: RMDir <directory path>

Arguments

<directory path>

This is a quoted string that defines the directory. Input the directory using the form <drive>/<dir></directory name>.<drive> and one or more <dir>s are optional.</directory name> stands for a directory name of up to 128 characters. While filename extensions are not required, they are highly recommended.

The current directory refers to the name of a directory as returned to the FILESystem:CWD? query.

Example

FILESystem:RMDir C:/"NOT-MINE"

This command creates the directory named NOT-MINE on the C drive.

FILESystem:WRITEFile

Description

This command (no query form) copies the GPIB port block data to a named file.

Group

File System

Syntax

FILESystem: WRITEFile <file path>, <data>

Arguments

• <file path>

This is a quoted string that defines the file name and path. Input the file path using the form <drive>/<dir>/<filename>.<drive> and one or more <dir>s are optional. <filename> stands for a filename of up to 128 characters. While filename extensions are not required, they are highly recommended.

• <data>

This is sent in the indefinite format of ARBITRARY BLOCK PROGRAM DATA encoding as described in IEEE 488.2.

Example

FILESystem: WRITEFile C:/"ABC.DEF", #0123<NL+EOL>

This command creates a file named "ABC.DEF" on the C drive, containing the three bytes (in hex format) 31, 32, and 33.